



## 5500sc Ammonia Monochloramine Standard 2 (2 ppm NH<sub>3</sub>), 2 L

**Product #:** 25238000  
**ZAR Price:** Contact Hach  
**Ships within 3 days**

Standard 2 for 5500sc AMC Analyzer, 2 L

---

### Specifications

Instrument:	5500sc
Measuring range:	N/A
Parameter:	Ammonia/Monochloramine
Storage conditions:	Keep container tightly closed in a dry and well-ventilated place. Keep in properly labeled containers.

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)