



Reagent 3 Package for 5500sc Silica Analyser

Product #: 6775102
ZAR Price: Contact Hach
Ships within 2 weeks

Includes 2 L bottle of deionized water and powder to be mixed together.

Specifications

Instrument:	5500sc Silica
Measuring range:	N/A
Parameter:	Silica
Storage conditions:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)