



5500sc Phosphate Low Range Reagent Set

Product #:

2035400

ZAR Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Complete set of reagents and standards for the 5500sc Low Range Phosphate analyser

Specifications

Instrument:	5500sc Phosphate
Measuring range:	N/A
Parameter:	Phosphate
Storage conditions:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)