



Be Right™



5500sc Standard 1 Low Range Phosphate 2L

Product #:

6776002

ZAR Price:

Contact Hach

2L bottle of Standard 1 for the 5500sc Low Range Phosphate analyzer

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)