



## 9245/9240 Sodium Standard, 100ppm, 1L

**Product #:** 2834253  
**ZAR Price:** Contact Hach  
**Ships within 3 days**

---

### Specifications

Concentration:	100 mg/L
Instrument:	9245 & 9240
Parameter:	Sodium
Quantity:	1 L
Storage conditions:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)