



## Sodium Standard Solution, 10 mg/L Na, 1 L

|                   |                |
|-------------------|----------------|
| <b>Product #:</b> | <b>2835153</b> |
| ZAR Price:        | Contact Hach   |
| Available         |                |

NaCl Standard, 10 mg/L Na, 1 L. For Polymetron 9240 & 9245 analyser.

---

### Specifications

|                     |                                                                                                 |
|---------------------|-------------------------------------------------------------------------------------------------|
| Concentration:      | 10 mg/L                                                                                         |
| Instrument:         | 9245 & 9240                                                                                     |
| Parameter:          | Sodium                                                                                          |
| Platform:           | Solution                                                                                        |
| Quantity:           | 1 L 1 L                                                                                         |
| Storage conditions: | Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture. |

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)