



## Sodium Nitrate, 0.5M, 500 mL

<b>Product #:</b>	<b>2507149</b>
ZAR Price:	Contact Hach
Ships within 1 week	

Sodium Nitrate, 0.5 M, 500 mL, for 9245 and 9240 Sodium analysers

---

### Specifications

Concentration:	5 mg/L
Instrument:	9245 & 9240
Parameter:	Sodium
Quantity:	500 mL
Storage conditions:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)