



## Hardness Indicator Solution, 5 - 100 mg/L, 500 mL

**Product #:** 2769249  
**ZAR Price:** Contact Hach  
**Ships within 3 days**

For use with the SP510 Hardness analyser.

---

### Specifications

Instrument:	SP510
Measuring range:	5 - 100 mg/L
Parameter:	Hardness
Storage conditions:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)