



## Magnesium Sulfate Standard, 29 mL, for Calcium Titration and Hardness Monitor, 10,000 mg/L as $\text{CaCO}_3$ , for SP510

**Product #:** 102233  
**ZAR Price:** Contact Hach  
**Ships within 3 days**

**For checking Hardness Monitor response: 10,000  $\pm$ 1000 mg/L as  $\text{CaCO}_3$**

---

### Specifications

Concentration:	10,000 mg/L
Instrument:	SP510
Parameter:	Hardness
Storage conditions:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)