



## Chromic Acid Cleaning Solution 500 mL for cleaning glassware.

Product #: 123349

 **Hazardous**

**OBSOLETE ITEM**  
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

Suggested replacements:  
• 203832

If necessary, we will change your selected shipping method to accomodate these items.

Chromic Acid Cleaning Solution 500 mL for cleaning glassware.

---

### Specifications

Instrument:	SP510
Measuring range:	N/A
Parameter:	Hardness
Storage conditions:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)