



APA6000 Hardness Standard Solution, 0.50 mg/L CaCO₃ (NIST), 946 mL

Product #:

2058016

ZAR Price:

Contact Hach

Ships within 3 days

Specifications

| | |
|---------------------|---|
| Concentration: | 0.50 mg/L |
| Description: | Hardness Standard Solution (NIST) |
| Instrument: | APA6000 |
| Parameter: | Hardness |
| Storage conditions: | Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture. |
| Volume: | 946 mL |

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)