



## CA610 Fluoride Standard 1, 0.5 mg/L, 473 mL

**Product #:** 2743811  
**ZAR Price:** Contact Hach  
**Ships within 3 days**

---

### Specifications

|                     |                                                                                                 |
|---------------------|-------------------------------------------------------------------------------------------------|
| Concentration:      | 0.5 mg/L                                                                                        |
| Instrument:         | CA610                                                                                           |
| Parameter:          | Fluoride                                                                                        |
| Storage conditions: | Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture. |

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)