



CA610 Fluoride Standard 2, 5.0 mg/L, 473 mL

Product #: 2797111
ZAR Price: Contact Hach
Ships within 3 days

Specifications

Concentration:	5.0 mg/L
Instrument:	CA610
Parameter:	Fluoride
Storage conditions:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)