



## CA610 Fluoride Reagent 1 (TISAB), 473 mL

**Product #:** 2812811  
**ZAR Price:** Contact Hach  
**Ships within 3 days**

---

### Specifications

Instrument:	CA610
Measuring range:	N/A
Parameter:	Fluoride
Storage conditions:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)