



CA610 Fluoride Reagent Set

Product #:	2816900
ZAR Price:	Contact Hach
No availability shown	

CA610 Fluoride Analyser Reagent Set includes: Reagent 1 (TISAB), Standard #1, 0.5 mg/L F, Standard #2, 5.0 mg/L F.

Specifications

Instrument:	CA610
Measuring range:	N/A
Parameter:	Fluoride
Storage conditions:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)