



Sedimeter for AISE, NISE and AN-ISE sc

Product #: LZY771
ZAR Price: Contact Hach

The sedimenter can be used as a bypass solution. The sample can be let through the sedimenter without prefiltering.

Specifications

Instrument:	suitable for:
	all AISE sensors;
	all NISE sensors;
	all AN-ISE sensors.
Length:	Intallation plate for wall mount: 800 mm x 495 mm
Manual Languages:	Non verbal instruction sheet
Material:	PVC
Mounting:	Wall installation (indoor only).
Weight:	Sedimeter incl. wall plate only: 14.7 kg