



Amtax Expulsion Solution 237 mL

Product #:

2825431

ZAR Price:
Available

Contact Hach



Hazardous

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Amtax Expulsion Solution 237 mL. The expulsion solution is used during the analyzer's automatic cleaning cycle. The cycle frequency is programmable using the Settings menu and can be manually activated from the Service menu. The Amtax Compact Analyzer requires three reagents: an expulsion solution, a standard solution and an indicator solution. Order a complete reagent set, which contains all three, based on parameter range: 2830700 (0.2-12 mg/L), 2830800 (2-120 mg/L) or 2830900 (20-1200 mg/L).

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)