



Sedimeter for Nitratax sc and Uvas sc

Product #:	LZX450
ZAR Price:	Contact Hach

The sedimenter can be used as a bypass solution. The sample can be let through the sedimenter without prefiltering.

Specifications

Instrument:	Suitable for all Nitratax and Uvas sc sensors
Length:	Intallation plate for wall mount: 800 mm x 495 mm
Manual Languages:	Non verbal instruction sheet
Material:	PVC
Pressure Range:	Not pressurized, open drain
Weight:	Sedimeter incl. wall plate only: 14.7 kg