



Reagent Set, Phosphax Compact Phosphate Analyzer

Product #:

2831000

ZAR Price:

Contact Hach

No availability shown



Hazardous

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Replacement reagent set for Hach Phosphax Compact Phosphate Analyzer.

Specifications

Range: 0.0 - 50.0 mg/L

Special Features: Phosphax Reagent Set

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)