



Optiquant SST Calibration Kit

Product #:	5733000
ZAR Price:	Contact Hach
Ships within 1 week	

Calibration kit, OptiQuant T-Line. Includes Calibration Cylinder (5737600), two 500 mL bottles of 800 NTU StablCal Standard (2660549), and probe bracket for Solitax sc immersion probes.

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)