



ORP buffer solution, 200 mV, 500 mL

Product #: 25M2A1001-115
ZAR Price: Contact Hach
Ships within 3 days

ORP Reference Solution used for ORP sensor calibration (in resealable plastic bottles).

Specifications

Instrument:	Digital Differential ORP
Measuring range:	200 mV
Parameter:	ORP
Quantity:	500 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)