



ORP buffer solution, 600 mV, 1L

Product #: 25M2A1002-119
ZAR Price: Contact Hach
Ships within 3-5 weeks

ORP Reference Solution used for ORP sensor calibration (in resealable plastic bottles).

Specifications

Instrument:	Digital Differential ORP
Measuring range:	600 mV
Parameter:	ORP
Quantity:	1.0 L

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)