



ORP buffer solution, 600 mV, 3.5 L

Product #:	25M2A1002-123
ZAR Price:	Contact Hach
Available	

ORP Reference Solution used for ORP sensor calibration (in resealable plastic bottle).

Specifications

Instrument:	Digital Differential ORP
Measuring range:	600 mV
Parameter:	ORP
Quantity:	3.5 L

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)