



Standard Cell Solution, Concentrated pH 7.0 Buffer (Equi-Transferrant), 500 mL

Product #: 25M1A1025-115
ZAR Price: Contact Hach
Ships within 3 days

Standard Cell Solution for Hach pHD Differential pH and ORP Sensors, packaged in resealable 500 mL bottle. Used to replenish standard cell chamber while replacing salt bridge.

Specifications

Concentration:	pH 7.0
Instrument:	pHD
Parameter:	pH
pH Range:	7 buffer solution
Quantity:	500 mL
Volume:	500 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)