



Be Right™



Refill kit for 32313

Product #: 32314

OBSOLETE ITEM
This item is no longer available.

Reagents to prepare a reference solution to calibrate CO₂ analysers. High carbonation 3 - 8 g/kg. Refill kit for 32313.

Traceable and Trustworthy Results

Effective calibration is critical to the performance of any analytical instrument and the ORBISPHERE 3231X is designed to give a traceable and trustworthy calibration standard with defined carbon dioxide content to $\pm 0.5\%$. Each reagent kit contains a certificate of analysis.

Quick and Easy to Use

Specifications

Instrument:	Orbisphere 3658 Analyzer Portable CO ₂
Measuring range:	3 - 8 g/kg
Parameter:	CO ₂
Storage conditions:	Store in a dry place at room temperature

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)