



Refill kit for 32315

Product #: 32316

OBSOLETE ITEM
This item is no longer available.

Reagents to prepare a reference solution to calibrate CO₂ analysers. Low carbonation 0 - 3 g/kg.

Specifications

Instrument:	Orbisphere 3658 Analyzer Portable CO2
Measuring range:	0 - 3 g/kg
Parameter:	CO2
Storage conditions:	Store in a dry place at room temperature

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)