



Recharge kit of 4 pre-filled cartridges with premounted 2952A membranes for A1100, GA2400 & GA2800 EX oxygen sensors

Product #:

2952A-A

ZAR Price:
Available

Contact Hach



Hazardous

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Includes O-rings, cotton filters and cleaning tools.

Specifications

Application:	Corrosion control; in-line beverage; deaerated water
Compatible Instruments:	A1100, GA2400 & GA2800 EX oxygen sensors
Material:	Membrane: ETFE
Membrane:	2952A membranes