



## Recharge kit of 4 pre-filled cartridges with premounted 29552A membranes for A1100, GA2400 & GA2800 EX oxygen sensors

Product #:

29552A-A

ZAR Price:

Contact Hach

No availability shown



**Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Includes O-rings, cotton filters and cleaning tools.

---

### Specifications

Application:

In-line wort; air/O<sub>2</sub> injection; sewage treatment

Material:

Membrane: PTFE