



Recharge kit of 4 pre-filled cartridges with premounted 2956A membranes for A1100, GA2400 & GA2800 EX oxygen sensors

Product #:

2956A-A

ZAR Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Includes o-rings, cotton filters and cleaning tools.

Specifications

Application:

Corrosion control; in-line beverage; deaerated water

Material:

Membrane: PFA