



Recharge kit of 4 pre-filled cartridges with premounted 2935A membranes for A1100, GA2400 & GA2800 EX oxygen sensors

Product #:

2935A-A

ZAR Price:

Contact Hach



Hazardous

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Includes o-rings, cotton filters and cleaning tools.

Specifications

Compatible Instruments:

A1100, GA2400 & GA2800 EX oxygen sensors

Membrane:

2935A membranes