



Recharge kit of 4 premounted 2956A membranes for A1100 oxygen sensors with 33051-xT cap only

Product #:

2956A-AT

ZAR Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Includes electrolyte, o-rings, cotton filters and cleaning tools.

Specifications

Instrument:	Orbisphere A1100-S00T O2 sensor, Orbisphere GA2x00-S00T O2 sensor
Parameter:	O2
Storage conditions:	Store in a dry place at room temperature

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)