



## Recharge kit of 4 pre-filled cartridges with premounted 2956A membranes for A1100 oxygen sensors with 33051-xT cap only

Product #:

2956A-AT05

ZAR Price:

Contact Hach



**Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Includes o-rings and cleaning tools. Special low concentration electrolyte.

---

### Specifications

Instrument:	Orbisphere A1100-S00T O2 sensor, Orbisphere GA2x00-S00T O2 sensor
Parameter:	O2
Storage conditions:	Store in a dry place at room temperature

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)