



Maintenance Kit for H₂ thermal conductivity sensors / membrane 29561A

Product #: 32739
ZAR Price: Contact Hach

Maintenance Kit for high pressure hydrogen thermal conductivity sensor. Includes membranes 29561A, membrane holding ring 29228.01, kit 29082, and tools for sensor maintenance.

Specifications

Instrument:	Orbisphere 31xxx TCD sensor
Parameter:	H2S
Storage conditions:	Store in a dry place at room temperature

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)