



Maintenance Kit for carbon dioxide, nitrogen, & hydrogen thermal conductivity sensors for use with the 3654 and 3658

Product #: 32746
ZAR Price: Contact Hach
No availability shown

Includes membranes 29561A, membrane holding ring 29228.01, O-rings 29039.0 and 28613.0, and tools for sensor maintenance

Specifications

Instrument:	Orbisphere 3654 Portable N2 and H2 Analyzer (TC), Orbisphere 3658 Analyzer Portable CO2 Analyzer (TC)
Parameter:	CO2, N2, H2
Storage conditions:	Store in a dry place at room temperature

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)