



Orbisphere 3658 Portable CO₂ (TC) analyser with N₂ purge

Product #:

3658/418



Hazardous

OBSOLETE ITEM

This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

The Orbisphere 3658 portable system provides accurate, selective CO₂ measurement. This instrument is specifically designed for the beverage market and is used in bottled water and soft drink plants, breweries and wineries.

Accurate measurements for low carbonated and very fizzy drinks and still wines.

Accuracy of ±1% of reading within ±5 °C of calibration temperature

Short analysis time, ideal for measurement in small packages

Response time of 1 minute to 90% of complete measurement.

Low maintenance - annual membrane change and calibration

Specifications

Parameter:	CO ₂
What's included?:	Orbisphere 3658 Portable CO ₂ Analyser, Thermal Conductivity Sensor, Communications Software, Two C-Type Batteries, User Manual

What's included?

Orbisphere 3658 Portable CO₂ Analyser, Thermal Conductivity Sensor, Communications Software, Two C-Type Batteries, User Manual