



Recharge kit of 4 pre-filled cartridges with premounted 29552A membranes for C1100 ozone sensors with 33051-xT cap only

Product #:

29552A-CT

ZAR Price:

Contact Hach

No availability shown

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Includes o-rings and cleaning tools.

Specifications

Instrument:	Orbisphere C1100 O3 sensor
Parameter:	O3
Storage conditions:	Store in a dry place at room temperature