



Maintenance kit for carbon dioxide & nitrogen TC sensors with protection cap 29111

Product #: 32760
ZAR Price: Contact Hach
No availability shown

Includes membranes 29561A, membrane holding ring 29228.01, 28002 washers, and tools for sensor maintenance.

Specifications

Instrument:	Orbisphere 31xxx TC Sensor
Parameter:	CO ₂ , N ₂
Storage conditions:	Store in a dry place at room temperature

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)