



Be Right™



Cylinder, Graduated 5 mL

Product #: 50897
ZAR Price: Contact Hach

- Class A Cylinders with Hex Base, Unserialised
- Double Metric Scale,
- Calibrated to contain (TC)
- Made to meet ASTM E-1272 norms
- Printed in durable white enamel

Specifications

| | |
|---------------------|---------------|
| Body Material: | Glass |
| Capacity: | 5 mL |
| Certifications: | Class A |
| Measurement method: | To Contain |
| Quantity: | 1 each |
| Type: | Double Metric |