



Flow Through Sample Chamber with Trap, 1” Pipe, for BioTector B7000i/B7000i Dairy

Product #: 19-BAS-014

ZAR Price: Contact Hach

Developed for samples containing sand, grit or heavy particulates, the Flow Through Sample Chamber with Sand Trap is a self-cleaning, simple device which does not require filters or screens. This sampler is designed to operate in conjunction with the BioTector sample pump. The body of the sampler is built from non-transparent material so as to prevent algae growth. The sample head is manufactured using transparent material for easy inspection. Algae growth on the sample head is prevented by a non-transparent sliding cover. The sample chamber is sealed and cannot contribute to build up of gas or vapours in a confined space.