



Intellical PHC805 Laboratory General Purpose Refillable Glass pH Electrode with Calibration and Maintenance Reagents Pack

Product #:
ZAR Price:

PHC80501AP
Contact Hach



Hazardous

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

The RedRod Sensors Now Bundled with Calibration Reagents

The RedRod sensor kit comes with calibration buffers reagents and maintenance solutions to provide an all-in-one easy to use reliable solution for your calibrations, helping you ensure accurate pH readings with your HQD pH Meter. The 500 mL bottles allows up to 6 months of pH electrode calibration when calibrating once per week, or up to 25 calibrations at a very attractive price compared if purchased separately.

Benefit from a significantly reduced price on calibration reagents

RedRod Sensors purchased as a kit together with reagents offer a 75% discount on reagents compared if purchased separately.

Reliable calibrations for optimum sensor performance

Bundle includes NIST traceable colour-coded pH 4.01, pH 7.00 and pH 10.01 buffer solutions and electrode cleaning solution for regular maintenance and pH electrode storage solution.

RedRod technology provides fast, accurate and reproducible measurements

Symmetry between internal and external RedRod reference elements ensures that the iso pH is the same as the zero pH, providing highly reproducible results and a fast, accurate response, even in samples with fluctuating temperatures. RedRod technology also provides a stable signal by protecting the reference elements from light interference, and the saturated KCl salt-bridge aids in creating a stable liquid junction potential.

RedRod reliability when measuring challenging samples over wider temperature ranges

RedRod encapsulation of reference elements means no silver ions and less risk of clogged junctions. It also enables accurate measurements over a broader temperature range than standard electrodes.

Intellical digital probes can be moved between meters without the need to re-calibrate or re-enter measurement settings

Ideal for a multi-user environment with multiple HQD series laboratory and portable meters.

Specifications

Accuracy:	±0.02 pH
Application:	High Performance (0 to 80 °C)
Cable Length:	1 m
Connector:	HQD specific
Electrode Type:	Standard
Filling Solution:	3M KCl
Kit?:	Yes
Length:	200 mm

Material:	Sensor Body: Glass
Method Type:	Laboratory: RedRod / Refillable Reference Element
Operating temperature range:	0 - 80 °C
Parameter:	pH
Probe Junction :	Porous Pin
Product Kit:	Model: PHC805
	Accessories Included: Calibration and Maintenance Reagents Pack
Range:	0 - 14 pH
Reference system:	RedRod - Ag/AgCl
Sensor material:	Glass
Special Features:	RedRod Technology Reference Element
Temperature Range:	0 to 80 °C
Test requirements:	Parameter Needed: pH
	Minimum Sample Depth (mm): 18
Warranty:	12 months
Water Resistance:	IP65
What's included?:	Intellical PHC805 Laboratory pH electrode with storage soaker bottle, 1 m cable, 5-mL bottle reference electrolyte filling solution KCl 3M solution), calibration certificate, and basic user manual. Buffer solution kit, colour-coded, pH 4.01, pH 7.00 and pH 10.01, 500 mL (#2947600). Electrode cleaning solution for regular maintenance, 500 mL (#2965249). pH electrode storage solution, 500 mL (#2756549).

What's included?

Intellical PHC805 Laboratory pH electrode with storage soaker bottle, 1 m cable, 5-mL bottle reference electrolyte filling solution KCl 3M solution), calibration certificate, and basic user manual. Buffer solution kit, colour-coded, pH 4.01, pH 7.00 and pH 10.01, 500 mL (#2947600). Electrode cleaning solution for regular maintenance, 500 mL (#2965249). pH electrode storage solution, 500 mL (#2756549).