



AT1000 Syringe Volume Verification Kit

Product #:	LZE228
ZAR Price:	Contact Hach
No availability shown	

Syringe volume verification kit for Titralab AT1000 Series automatic titrator

Specifications

Warranty:	12 months
What's included?:	Kit includes: stand, boss head and clamp (LZE231); tubing, pipet tip, pipet holder and tube-syringe holder (LZE230); syringe verification activation license (LZE229)

What's included?

Kit includes: stand, boss head and clamp (LZE231); tubing, pipet tip, pipet holder and tube-syringe holder (LZE230); syringe verification activation license (LZE229)