



CL17sc Cell Cleaning Kit

Product #: 8573100
ZAR Price: Contact Hach
No availability shown



Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Hach recommends cleaning the CL17sc cell on a monthly (or as needed) basis. Frequency of cleaning may depend on the application and sample qualities.

The CL17sc Cell Cleaning Kit contains 5.25N sulfuric acid (100 mL dropper bottle) and cotton swabs (10).

Specifications

Instrument:	CL17sc
Measuring range:	N/A
Parameter:	N/A
Storage conditions:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.
What's included?:	Sulfuric acid (100 mL dropper bottle) and cotton swabs (10).

What's included?

Sulfuric acid (100 mL dropper bottle) and cotton swabs (10).