



Ammonia TNTplus Vial Test, HR+ (10 - 100 mg/L $\text{NH}_3\text{-N}$), 25 Tests

Product #:

TNT837

ZAR Price:

Contact Hach



Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Expert High Range + Ammonia test made simple.

For samples with high ammonia concentrations - ideal for monitoring ammonia concentrations of wastewater inlet samples.

Ammonia TNTplus chemistry and your Hach spectrophotometer are engineered to simplify water analysis for accurate results, everytime. For determination of high concentrations of ammonia, e.g. in municipal or industrial wastewater inlets by the salicylate method.

With Ammonia TNTplus reagents, simply add sample and add the Dosicap zip reagent, shake the vial and insert into your Hach spectrophotometer! The Direct Read Technology will provide your ammonia results in mg/L! No more sample dilutions, cleaning glassware, dealing with large amounts of waste, and worrying about human errors.

25 Ammonia vials per box.

Easy and safe handling

The innovative vial and reagent delivery makes tests easier to use than powder pillows or liquid reagents and significantly reduces spillage or contamination risk.

No reagent blank necessary

The high quality of TNTplus vials, tight reagent production controls, instrument calibration verification, and high instrument stability all combine to eliminate the need to run reagent blanks—saving you time and money!

Automatic method detection

The bar-coded vials allow for automatic method detection and measurement by the spectrophotometer, considerably increasing the speed of analysis. Errors are reduced as instrument averages 10 readings and rejects outliers from scratched, flawed, or dirty glassware.

Documented shelf life and COA

The 2D vial barcode details batch number and expiry date of reagents, which are documented along with the measurement result. An automatic warning is issued if expiry date has passed. Lot specific information (COA) is available on RFID tag on the box.

Specifications

Instrument:	DR3900, DR6000, DR1900
Method:	10205
Number of tests:	25
Parameter:	Ammonia, Nitrogen
Platform:	TNTplus
Range:	10 - 100 mg/L $\text{NH}_3\text{-N}$
Storage conditions:	2 - 8 °C (keep refrigerated)

Truecal:

Yes