



## Ammonium Standard Solution 50 mg/L $\text{NH}_4\text{-N}$ , 100 mL

**Product #:** LCA503  
**ZAR Price:** Contact Hach

Ammonium Standard Solution with 50 mg/L  $\text{NH}_4\text{-N}$  for quality control checks with LCK503 Ammonium cuvette test. LCA503 is also suitable for quality checks with the Amtax sc Ammonium analyser.

---

### Specifications

Concentration:	50 mg/L $\text{NH}_4\text{-N}$
Parameter:	Ammonium
Quantity:	100 mL
Storage conditions:	15 - 25 °C (protect from light)