



Ammonium Standard Solution 2.5 mg/L $\text{NH}_4\text{-N}$, 100 mL

Product #: LCA505
ZAR Price: Contact Hach

Ammonium Standard Solution with 2.5 mg/L $\text{NH}_4\text{-N}$ for quality control checks with LCK505 Ammonium cuvette test. LCA505 is also suitable for quality checks with the Amtax sc Ammonium analyser.

Specifications

Concentration:	2.5 mg/L $\text{NH}_4\text{-N}$
Parameter:	Ammonium
Quantity:	100 mL
Storage conditions:	15 - 25 °C (protect from light)