



## USB + Power Module for 2100Q/2100Q IS & DR 1900

**Product #:** LZV813.99.00002

**ZAR Price:** Contact Hach  
**Available**

The USB + Power Module is an optional accessory for the 2100Q & 2100Q IS portable turbidimeters, and the portable spectrophotometer DR1900. This module adds further functionality to the instrument allowing the user to:

**Transfer data to a PC in XML format**

**Transfer data to the Citizen PD-24 printer (2960100)**

**Operate the instrument on line power (100-240Vac, 50/60 HZ) and Charge NiMH batteries (2971304)**

**Perform firmware updates**

**Quickly enter Sample and User IDs with the optional barcode scanner (LZV566)**

---

### Specifications

Includes:	Includes module, module cover, power supply with four plug adapters, standard USB cable with Mini-USB connector, special USB cable adapter, and manual.
Input Voltage:	100 - 240 VAC, 50/60 Hz
Output Voltage:	9V
Warranty:	12 months