



## HR Manganese Chemkey Reagents (box of 25)

<b>Product #:</b>	<b>3007000</b>
ZAR Price:	Contact Hach
No availability shown	

### High Range Manganese Chemkey for SL1000 Portable Parallel Analyzer

For determination of Manganese using Hach Method 10319

Ensure SL1000 has been updated with Software Version V3.0.12.1489 to use this Chemkey.

#### Less Variability

Automation and internal temperature control make the entire process consistent and repeatable, while applying the same processes as current Hach methods. Avoid manual steps that can introduce variability, even when performed by experienced testers.

#### Less Hassle

All chemicals and processes are entirely contained inside the Chemkey. There are no powder pillows or glass vials to handle.

---

### Specifications

Instrument:	SL1000
Number of tests:	25
Parameter:	Manganese, HR
Quantity:	pk/25
Range:	0.10 - 25.0 mg/L Mn