



Reagent Kit for EZ1010 Copper Cu(II) Analyzer

Product #:

APPC1010KTO



Hazardous

ZAR Price:

Contact Hach

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Specifications

Parameter:

Copper, dissolved

What's included?:

APPC101001 Nitric Acid Solution, 0.25 M, 2 L (2x)

APPC230201 EZ2302, Color Indicator Solution, 2 L (2x)

259342 Standard Solution, Copper, 1000 mg/L Cu (NIST), 100 mL (1x)

What's included?

APPC101001 Nitric Acid Solution, 0.25 M, 2 L (2x)APPC230201 EZ2302, Color Indicator Solution, 2 L (2x)259342 Standard Solution, Copper, 1000 mg/L Cu (NIST), 100 mL (1x)