



Cylinder, Graduated, Mixing, Glass, 50 mL $\pm$ 0.4 mL, 1.0 mL divisions, glass stopper #16

Product #: 104592
ZAR Price: Contact Hach

Tall form, single metric scale, glass stopper, graduated. Calibrated "to contain," durable white scale.

Scale is durable white ceramic enamel. Designed to meet ASTM Specification E1272, Style 2 (Class B). Includes glass stopper.

Specifications

Body material:	Glass
Capacity:	50 mL
Certifications:	Class B
Measurement Method:	To Contain
Quantity:	1 each
Type:	Single Metric, Mixing