



Reagent Kit for EZ1028 Nitrite Analyser (Reagent 1 & 2)

Product #:

APPC1028KTO

ZAR Price:

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Ready-to-use Reagents for EZ Series Online Analysers

Ready-to-use

Hach EZ reagents are ready-to use with EZ Series Online Analysers. Since there is no additional preparation needed this eliminates an additional error source when using the EZ instruments and reagents system.

Confidence in your results

Reduce the sources of errors in your workflow to increase the confidence in your results.

Less complexity

Using ready-to-use EZ reagent kits from Hach can reduce the number of handling steps around your analytical equipment significantly which allows less hands-on-time and complexity.

Specifications

Parameter:	Nitrite
What's included?:	52022 Sulfanilamide Reagent 1, 400mL (10x)
	52023 Sulfanilamide Reagent 2 (10x)

What's included?

52022 Sulfanilamide Reagent 1, 400mL (10x) 52023 Sulfanilamide Reagent 2 (10x)